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

HP Fff – firmness of fruit pulp

List of Contents

List of Contents	3
1 General information	4
1.1 Safety precautions	4
1.2 Volume of delivery	4
2 Hardness Tester	5
2.1 Device Description	5
2.2 Ranges of Application	5
2.3 Technical Details	5
3 Installation	6
3.1 Control of contents	6
3.2 Installing and replacing the test anvil	6
4 Measurement	7
4.1 Measuring principle	7
4.2 Determination of fruit flesh firmness	7
4.3 Determination of flesh firmness on soft and small fruits	8
5 Servicing the measuring instrument	9
5.1 Conditions of warranty	9
5.2 Service and maintenance	9
5.3 Care	9
6 Procedure for troubleshooting	10
7 Recycling, Information and Contact Information	11
7.1 Recycling and Waste Disposal	11
7.2 Information for return of goods	11
7.3 Contact Information	11
8 Technical Data and Accessories	12
8.1 Technical Data	12
8.2 Accessories	12
9 EU-Declaration of conformity	13
Table of Figures	14
Company overview	15

1 General information

Although the information contained in these operating instructions was controlled carefully for accuracy and completeness, no liability can be taken for errors or omissions. These operating instructions may not be multiplied partly or completely in any kind or translated to another language without our previous written consent.
Keep for future application! Technical changes without notice!

1.1 Safety precautions

When working with the HP Fff, named in the following hardness tester, the following instructions must be observed:

Attention!

Works on hardness tester may only be done by authorized persons.

- The hardness tester may only be used for hardness determination on material as described under Ranges of Application, *see chapter 2.2, „Ranges of Application“, Page 5.*
- The hardness tester is to be sheltered from dusty, oily, greasy and metal-dusty air, sources of heating (direct sun beaming, ovens), humidity, wetness and vibration as well as from damage caused by falling down.
- A soft and lint-free cloth should be used for cleaning.

Attention!

Alcohol, gasoline, diluents or other easily inflammatory substances may not be used. The use of such substances can lead to fires.
The cleaning cloth should be soft and lint free.

1.2 Volume of delivery

See delivery note

2 Hardness Ttester

2.1 Device Description

- (1) Bottom side of the hardness tester
- (2) Pressure plate
- (3) Test anvil

FIG. 1 FRONT VIEW

2.2 Ranges of Application

Test anvil	Ranges of Application
Cylinder 0,10 cm2	Peaches, apricots, plums, damsons
Cylinder 0,25 cm2	cherries, melons, tomatoes, blueberries, grapes, aubergine
Cylinder 0,50 cm2	strawberries, watermelons, broccoli, meat, tofu
ball Ø 2,5 mm	cucumber, radish, carrot, kohlrabi, courgette, onion (firm)
ball Ø 5,0 mm	apples, avocados, pears, cucumbers, oranges, peppers, papayas, melons, mangoes, citrus fruits, onion (soft)
Segment of a sphere r=6,3 mm	cooked potatoes, bananas, smoked fish

2.3 Technical Details

Test method	Spring Force [cN]	Contact pressure force [N]	Test anvil	Measuring distance [mm]	Reading range
Fff - 0,1	806,50	12,5	0,10 cm²	2,5	0 - 100
Fff - 0,25			0,25 cm²		
Fff - 0,5			0,50 cm²		
Fff - 5,0			Ø 5,0 mm		
Fff - 2,5			Ø 2,5 mm		
Fff - r=6,35			Ø 12,7 mm		

English

3 Installation

3.1 Control of contents

Check the supplied equipment for completeness and soundness, see delivery note.

3.2 Installing and replacing the test anvil

Test anvil	drilling in the knurled screw
Cylinder 0,10 cm ² (Ø 3,65 mm)	Ø 4,1 mm
Cylinder 0,25 cm ² (Ø 5,64 mm)	Ø 6,1 mm
Cylinder 0,50 cm ² (Ø 7,98 mm)	Ø 8,5 mm
Ball Ø 2,5 mm	Ø 3,0 mm
Ball Ø 5,0 mm	Ø 5,5 mm
Spherical section r=6,3 mm	Ø 10,5 mm

Attention!

Please take care that you always use test anvil with fitting pressure plate. If not, you will have a faulty measurement.

- Select test anvil (3) and fitting pressure plate (2).
- Push the test anvil into the drilling of the pressure plate.

Attention!

The test anvil and the pressure plate must be cleaned regularly to avoid measuring errors because of gluing. Use isopropanol for this purpose.

Screw the pressure plate together with the test anvil into the thread (9) on the bottom side of the hardness tester.

Removing the test anvil, proceed in reverse order.


FIG. 2 TEST ANVIL AND PRESSURE PLATE


FIG. 3 FITTING THE TEST ANVIL

4 Measurement

4.1 Measuring principle

The yielding of the surface is indicated at a specified contact pressure of 1 kg, a defined measuring distance and a measuring time of 3 s.

The test anvil does not penetrate the fruit peel.

With the measuring device the deformation of the fruit is observed under a certain and constant contact pressure and a defined measuring distance.

The measured value indicates no force unit but a quotient.

Therefore, there is an absolute need to mention the test anvil which had been used for the hardness test

On each fruit 2 measurements are made on opposite sides.

The average value is calculated from the 2 measurements.

For denomination of the average value please use.:

e.g. middle index of firmness, 0.1 cm^2

4.2 Determination of fruit flesh firmness

- Make sure that the test anvil moves smoothly.

Attention!

The test anvils may only be used with the supplied hardness tester because they are matched to it.

- Press the hardness tester onto the fruit peel (4) and read the value after 3 seconds.

It is absolutely necessary to indicate the test anvil that was used for the hardness test.

Attention!

During the measurement, the pressure should not be applied too fast but smoothly, otherwise too high values will be displayed.

- Make a random test on 30 fruits.
- Make 2 measurements on the opposite sides on each fruit, e.g.: sun side and shadow side.

Attention!

The measuring point should be as smooth as possible and without skin defects.


FIG. 4 MEASUREMENT

4.3 Determination of flesh firmness on soft and small fruits

- Hold the measuring device vertically onto the soft or small fruit (5) so that the contact surface of the pressure plate rests on the fruit.

In this case, sufficient contact pressure is applied on the soft or small fruit.

For a reliable measurement the surface of the fruit skin or skin should be even.

This is especially true for soft and small fruits whose diameter is smaller than the diameter of the printing plate, such as strawberry, cherry, grape, etc.

Attention!

The fruit is destroyed when pressing down the measuring device.

FIG. 5 MEASUREMENT OF SOFT AND SMALL FRUITS

5 Servicing the measuring instrument

5.1 Conditions of warranty

Please see our common business conditions AGB's (www.bareiss.de)

Attention!

There is no claim of warranty for damages or faults caused by:

- ignoring the correct connection
- inappropriate handling
- neglecting the operating instructions
- repairing done by persons without authorization
- removing the type plates

5.2 Service and maintenance

A careful inspection and maintenance of the hardness tester is recommended.

When maintenance is performed by the user, the following operations are possible:

- If necessary, clean the hardness tester (see *chapter 5.3, „Care“, Page 9*).
- Contact the manufacturer if the hardness tester is damaged or has got some malfunctions (see *chapter 7.3, „Contact Information“, Page 11*).

5.3 Care

- Only mild cleaning agents may be used for cleaning the measuring device in order to avoid surface damages.
- The cleaning cloth must be soft and lint-free.

Attention!

Do not use alcohol, petrol, thinners or other highly inflammable substances.
The use of such substances can lead to fires.

6 Procedure for troubleshooting

Malfunction	Cause	Procedure
Measurement result incorrect	Damaged indenter Changed spring adjustment adjusted	Send device for repair
Hardness values are too high	Pressing too fast Irregular pressing	Press the hardness tester evenly
Hardness values deviate a lot	the test anvil is dirty	Clean the Test anvil -- see chapter 3.2, <i>Installing and replacing the test anvil</i> , Page 6

If the mentioned procedures do not correct the malfunction, a maintenance of the hardness tester is necessary.

For questions regarding service and troubleshooting,
please contact Bareiss Prüfgerätebau GmbH phone, fax or e-mail
(see chapter 7.3, „Contact Information“, Page 11).

7 Recycling, Information and Contact Information

7.1 Recycling and Waste Disposal

Please dispose of recyclable materials in an environmentally friendly manner. Old devices can be disposed of at suitable recycling collection points. If you cannot recycle the equipment yourself, the manufacturer of the equipment will do this for you. Send us your device with the note "Recycling device".

7.2 Information for return of goods

The hardness tester must be tested extensively before returning it, e.g. there could be a defect or malfunction due to a faulty application

(see chapter 6, „Procedure for troubleshooting“, Page 10).

In the case of uncertainties, you can contact us at any time by telephone / fax / e-mail.

To avoid further inquiries, a precise error description as well as the S/N are necessary.

For returns please use the "Service Form", which is available for download on our website.

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8 Technical Data and Accessories

8.1 Technical Data

Scaling	1
Test method	Fff
Measuring time	3 s
Dimensions in mm WxDxH	65 x 40 x 100
Weight	0,300 kg

8.2 Accessories

Denomination	Product no.
Accessories Functional tests	
Test anvil 0.10 cm ² with knurled screw	fm00015
Test anvil 0.25 cm ² with knurled screw	fm00016
Test anvil 0.50 cm ² with knurled screw	fm00017
Test anvil Ø 5.0 mm with knurled screw	fm00018
Test anvil Ø 2.5 mm with knurled screw	fm00019
Test anvil r = 6,35 mm with knurled screw	fm00021
Test stand BS 61 II	fm00029-1

9 EU-Declaration of conformity

EU – Declaration of conformity

according to the
EU-Measuring Instruments Directive
EC-Machinery Directive

2014/32/EU EU-Abl. L 96/149 of 26.02.2014
2006/42/EC EU-Abl. L 157/24 of 09.06.2006

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Authorized person for documentation: Mr. Harald Glöggler, see address of manufacturer

Product name:

Hardness tester

Type:

HP

Serialnr.:

see Type plate

We herewith confirm that the above-mentioned device has been developed, designed and manufactured in accordance with the guidelines, mentioned in this declaration.

Applied national standards and specifications:

DIN EN ISO 12100:2011

Machinery safety -
General principles for design -
Risk evaluation and Risk reduction

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2020-01-27

Date

Harald Glöggler
Authorized person for documentation

English

Table of Figures

FIG. 1 FRONT VIEW..... 5

FIG. 2 TEST ANVIL AND PRESSURE PLATE..... 6

FIG. 3 FITTING THE TEST ANVIL..... 6

FIG. 4 MEASUREMENT 7

FIG. 5 MEASUREMENT OF SOFT AND SMALL FRUITS 8

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Bareiss Shanghai Limited
Bareiss USA Inc.