

Ranges of Application / technical data

Test Method	Range of Application	Standards	Minimum material thickness [mm]
Abrasion Tester			
Abrasion resistance	The abrasion test serves for the evaluation of the resistance of elastomers against the rubbing wear. The wear condition is determined by the volume loss of a cylindrical specimen.	DIN 53516, DIN ISO 4649:2002	6*

* If no sample test specimens of the necessary thickness are available, the necessary specimen thickness can be achieved by gluing several specimens of minimum 2 mm thickness one upon the other. Their hardness should be at least 80 IRHD.

BAREISS