


THE WORLD OF
Material Testing
bareiss®

Overview Barofix & Centrofix

Centering devices for O-Rings, hoses, cables
and other finished products

Barofix Basic

- quick-center-device for O-rings
- for a cord diameter from $\varnothing$ 0,6 up to $\varnothing$ 7,0 mm
- outer diameter up to 50 mm
- Article No. fm11110



Barofix

- Quick centring device for O-rings
- for cord diameter from $\varnothing$ 0.6 to $\varnothing$ 8.0 mm
- outer diameter up to 100 mm
- Article No. fm11109

Optional hold-down
clamp for bigger O-Rings



Optional Plexiglas Plate available

Barofix 4

- automatic specimen centering device for method IRHD M, Micro Shore A)
- with optical specimen detection
- for O-rings from 1 mm to 10 mm cord thickness
- for plane-parallel plates
- thickness measurement
- 2-point measurement
- multiple measurement
- travel distance 90 mm
- support table Ø120 mm
- Article No. fm11106




Optional Plexiglas Plate available



Built-in thickness measurement function. Single-point or multi-point measurement, and multiple samples can be tested at the same time. Fully automatic positioning, *eliminating manual operation errors.*

Centrofix for hoses, cables and other parts


Article No. Centrofix Base fm01052



fm01052-04



fm01052-03



fm01052-02



Bareiss Prüfgerätebau GmbH
Germany

www.bareiss.de

