

# ROTOFIX

EN

Positioning device available as an accessory for the automatic measurement of specimens for the digi test II and the Gelomat. Different templates allow for the measurement of a varying number of samples and measuring positions.



The Rotofix is the smart automation solution for hardness testing on multiple samples with just one push of a button. In combination with the digi test II or the Gelomat, it provides the perfect entry point into automation.

# ROTOFIX

EN

## FEATURES

Through the electronic unit, the number and position of measurements can be freely defined. The measurement positions are set on the first sample and then automatically transferred to the other samples.

Thanks to our high vertical integration, it is possible to produce the template accurately fitting to your requirements.



Depending on the size of the samples, the number of samples on the template can be chosen in advance.

### TECHNICAL DATA

**Dim. W x H x D:** 26 x 14 x 10 cm

**Weight:** 2,82 kg

### SCOPE OF DELIVERY

**Rotofix**

**Connection cable electronic unit - Rotofix**

**Template**

# ROTOFIX

EN

## ACCESSORIES



**Templates for  
various sample sizes  
and geometries**

### REFERENCE

For the automatic testing of a larger number of samples,  
we recommend the Rotation XL.

MADE IN GERMANY SINCE 1954.

#### **Bareiss Prüfgerätebau GmbH**

DAKKS-Kalibrierlaboratorium

Breiteweg 1

89610 Oberdischingen, Germany

Tel +49 (0) 7305 / 96 42-0

Fax +49 (0) 7305 / 96 42-22

[sales@bareiss.de](mailto:sales@bareiss.de)



The accreditation is valid for the scope listed in  
certificate D-K-15206-01-00 (mechanical measurands  
in the range of hardness).