



# ROTOFIX FOR DIGI TEST II

EN

Optional Positioning Device for Automatic Measurement of Test Specimens with the digi test II. Different templates allow for the measurement of varying numbers of samples and measuring positions.



The Rotofix is the smart automation solution for hardness testing on multiple samples at the push of a button. In combination with the digi test II, it offers the perfect entry point into automation.

# ROTOFIX FOR DIGI TEST II

EN

## KEY FEATURES




### Flexible Sample Quantity

Depending on the size of the samples, the number of samples on the template can be pre-selected.



### Freely Definable Number of Measuring Points

The number of measuring points can be freely defined via the electronic control unit.



### Custom Sample Templates

Thanks to our high level of in-house manufacturing, we can produce templates precisely tailored to your requirements.



### Operation of the Rotofix via BareissOne Software

From creating a test plan to automated data acquisition during measurement, the software offers a wide range of possibilities.

## TECHNICAL SPECIFICATIONS

Measurements (W x D x H) 20 x 14 x 13 cm

Weight 3,99 kg

MADE IN GERMANY SINCE 1954.

**Bareiss Prüfgerätebau GmbH**  
DAKKS-Kalibrierlaboratorium  
Breiteweg 1  
89610 Oberdischingen, Germany  
Tel +49 (0) 7305 / 96 42-0  
Fax +49 (0) 7305 / 96 42-22  
sales@bareiss.de



The accreditation is valid for the scope listed in certificate D-K-15206-01-00 (mechanical measurands in the range of hardness).