

ROTOFIX FOR GELOMAT

EN

Optional positioning device for the Gelomat, enabling automatic measurement of gel capsules, gelatin-like materials, and culture media. The universal template allows for the measurement of varying numbers and sizes of samples.



The Rotofix is the smart automation solution for hardness testing on multiple samples at the push of a button. In combination with the Gelomat, it offers the perfect entry point into automation.

ROTOFIX FOR GELOMAT

EN

KEY FEATURES



Flexible Sample Quantity

The maximum number of samples is 10, but it can be flexibly adjusted to meet specific requirements.



Freely Definable Number of Measuring Points

The number of measuring points can be freely defined via the electronic control unit.




Reduced User Influence Thanks to Setting Gauges

The included setting gauges minimize user influence and ensure consistent operation of the sample template.



Custom Sample Templates

Thanks to our high level of in-house manufacturing, we can produce templates precisely tailored to your specific requirements.



Operation of the Rotofix via Celeste Software

From creating a test plan to automated data acquisition during measurement, the software offers a wide range of functionalities.

TECHNICAL SPECIFICATIONS

Measurements (W x D x H) 20 x 14 x 18 cm

Weight 4,87 kg

Bareiss Prüfgerätebau GmbH
DAKKS-Kalibrierlaboratorium
Breiteweg 1
89610 Oberdischingen, Germany
Tel +49 (0) 7305 / 96 42-0
Fax +49 (0) 7305 / 96 42-22
sales@bareiss.de



MADE IN GERMANY SINCE 1954.



The accreditation is valid for the scope listed in certificate D-K-15206-01-00 (mechanical measurands in the range of hardness).